


HLC-723G11 **B-THALASSAEMIA ANALYSIS MODE** THE ULTIMATE HAEMOGLOBINOPATHY SOLUTION FOR YOUR LAB



TOSOH BIOSCIENCE

Haemoglobinopathies are among the most common hereditary diseases of the world's population. About 4.5% of all human beings carry a gene for a thalassaemia or haemoglobin anomaly⁽¹⁾. The areas in which such abnormalities were originally most common extend from Africa over the Mediterranean basin and the Near- and Middle East to Southeast Asia and the Indian subcontinent. Global migration in the modern period has led to a continual spread of these anomalies to all regions of the world, with the result that they are rapidly becoming more common in the industrialised regions of Northern and Central Europe as well⁽²⁾.



Variants of thalassaemias and main abnormal haemoglobins interact to produce a wide range of clinical disorders of varying severity⁽³⁻⁴⁾. Homozygotes for ß-thalassaemia may develop either thalassaemia major or thalassaemia intermedia. Individuals with thalassaemia major are usually diagnosed within the first 2 years of life and require regular blood transfusions to survive⁽⁶⁾.

LABORATORY DIAGNOSIS

Diagnosis of Beta-Thalassaemia and other types of haemoglobinopathies should be done based on clinical symptoms (if available) and a number of laboratory tests, such as MCV, MCH, total red cell count, HbF, HbA2 and iron markers.


TOSOH G11 BETA-THALASSAEMIA SOLUTION







As a guideline, the below scheme can be used (Adapted from Mosca et at. J.Clin.Pathol. 2009 - with permission)



MLPA: Multiplex ligation-dependent probe amplification

- Quantitative determination of HbF and HbA, in 5 minutes • Chromatographic separation between HbA, and HbE
- High resolution chromatogram thanks to Tosoh's over
 - 40 years' experience in HPLC
- Full reagent traceability
- Easy to use and intuitive instrument
- Highly reliable system
- Instrument connectable to open laboratory automation lines

reased - 8.0%	
ormal or	
reased	
1-8%	
ormal	
h HbA ₂ al carrier	
*	
quence	

- Software to assist in the interpretation of chromatograms
- Tosoh Bioscience offers an on-line clinical interpretation platform under the supervision of highly knowledgeable experts in the field

Example chromatograms: Normal, &-Thalassaemia trait, HbS trait, HbE trait

** THALASSEMIA REPORT **	** THALASSEMIA REPORT **	** THALASSEMIA REPORT **	** THALASSEMIA REPORT **
G11 - 61461 2017/08/22 09:15 0P:ANN V03.00 N0: 0001 SL 0001 - 01 ID: 001052 CAL F Y= 0.9382X - 0.0044 A2 Y= 1.4332X + 0.2912	611 - 61461 2017/08/11 10:47 0P:ANN V03.00 N0: 0009 SL 0001 - 03 ID: 5373729 CAL F Y= 0.9362X - 0.0044 A2 Y= 1.4332X + 0.2912	G11 - 61461 2017/08/22 09:20 OP:ANN V03.00 N0: 0002 SL 0001 - 02 ID: 011052 CAL F Y= 0.9382X - 0.0044 A2 Y= 1.4332X + 0.2912	G11 - 61461 2017/06/27 12:00 0P:ROOT V03.00 N0: 0003 SL 0001 - 03 ID: 0001 - 03 CAL F Y= 0.9245X + 0.0306 A2 Y= 1.4358X + 0.3809
TP 829	TP 1320	TP 1028	TP 1615
Column: 2017/06/30 Buff.1: D1-101C 2017/08/22 Buff.2: D1-201C 2017/07/28 Buff.3: D1-301C 2017/07/27 H&W : HW065G/B 2017/07/26 CalSet: BCALO1C 2017/07/19	Column: 2017/06/30 Buff.1: D1-101C 2017/07/28 Buff.2: D1-201C 2017/07/28 Buff.3: D1-301C 2017/07/27 H&W : HW065G/B 2017/07/26 CalSet: 1BCALD1C 2017/07/19	Column: 2017/06/30 Buff.1: D1-101C 2017/08/22 Buff.2: D1-201C 2017/07/28 Buff.3: D1-301C 2017/07/27 H&W : HW0656/B 2017/07/26 CalSet: 1BCAL01C 2017/07/19	Column: 2017/06/26 Buff.1: 2017/06/28 Buff.2: 2017/06/26 Buff.3: 2017/06/28 Buff.3: 2017/06/29 L&W : 2017/06/29 CalSet: 2017/06/26
NAME % TIME AREA F 2.08 0.61 62.77 A0 77.06 1.81 2178.03 A2 2.17 2.16 37.14 E+ 0.00 0.00 0.00 D+ 0.00 0.00 0.00 S+ 0.00 0.00 0.00 C+ 0.00 0.00 0.00 TOTAL AREA 2826.50	NAME % TIME AREA F 0.71 0.55 13.38 A0 76.26 1.81 1347.57 A2 4.91 2.14 56.90 E+ 0.00 0.00 0.00 D+ 0.00 0.00 0.00 S+ 0.00 0.00 0.00 C+ 0.00 0.00 0.00 C+ 0.00 0.00 0.00 C+ 0.00 0.00 0.00 T0TAL AREA 1767.05	NAME % TIME AREA F 5.09 0.61 144.03 AO 44.02 1.82 1168.57 A2 5.71 2.18 101.35 E+ 0.00 0.00 0.00 D+ 0.00 0.00 0.00 S+ 26.29 3.07 697.78 C+ 0.00 0.00 0.00 TOTAL AREA 2654.46 46	NAME % TIME AREA F 0.71 0.55 7.08 AO 65.35 1.82 630.88 A2 3.23 2.14 20.48 E+ 11.66 2.67 112.54 D+ 0.00 0.00 0.00 S+ 0.00 0.00 0.00 C+ 0.00 0.00 0.00 C+ 0.00 0.00 0.00
F : 2.08% A2: 2.17%	F : 0.71% A2: 4.91%	F : 5.09% A2: 5.71%	F : D.71% A2: 3.23%
P00 2.30 0.21 64.95 P01 1.05 0.31 29.74 P02 0.12 0.74 3.25 P03 1.64 1.00 46.48 P04 5.26 1.21 148.56 P05 5.69 1.30 160.70 P06 1.70 2.69 48.18 P07 0.26 3.32 7.47 P08 1.39 3.77 39.23	$\begin{array}{cccccccccccccccccccccccccccccccccccc$	$\begin{array}{cccccccccccccccccccccccccccccccccccc$	$\begin{array}{cccccccccccccccccccccccccccccccccccc$

References

- 1. Angastiniotis M, Vives Corrons J-L, Soteriades ES, Eleftheriou A. The Impact of Migrations on the Health Services for Rare Diseases in Europe: The example of Haemoglobin Disorders, in: The Scientific World Journal Volume 2013
- 2. Haemoglobinopathies on the move: Is Europe ready? Report by group of experts from the European Network for Rare and Congenital Anaemias (ENERCA) and the Thalassemia International Federation (TIF) in collaboration with the International Organization for Migration (IOM), Migration Health Division, Regional Office Brussels

3. Weatherall DJ, Clegg JB. 2001. The thalassaemia syndromes. 4th Edition. Oxford: Blackwell Science Ltd 2001

4. Thein SL. The molecular basis of ß-thalassemia. Cold Spring Harb Perspect Med. 2013 May 1;3⁽⁶⁾

5. Cao A, Kan YW. The prevention of thalassemia. Cold Spring Harb Perspect Med. 2013 Feb 1;3⁽²⁾



 EC
 REP

 Transportstraat 4- 3980 Tessenderlo- BELGIUM

 Tel: +32 (0)13 66 88 30- Fax: +32 (0)13 66 47 49

sales-marketing.tsen@tosoh.com www.tosohbioscience.eu

"HLC" and "HLC-723" are the registered trademark of Tosoh Corporation in EU, etc.

