



Automated Immunoassay Analysers

# AIA<sup>®</sup>-360

## The Special One



Automated Immunoassay Analyser AIA<sup>®</sup>-360

Better diagnostics **for better living**

**TOSOH BIOSCIENCE**

# AIA<sup>®</sup>-360

## The Special One



<b>METHODOLOGY</b>	Reagent type	AIA-Pack <sup>®</sup>
	Test principle	1 AIA-Pack <sup>®</sup> = 1 test
	Immunochemical determination	Entirely inside AIA-Pack <sup>®</sup>
	Signal detection	Bichromatic Fluorescence Kinetics
<b>CAPACITY</b>	Time to first result	18 minutes
	Sample loading capacity	25
	Number of different tests per sample	4
	Test loading capacity	25 (if 1 per sample)
	Calibration stability	90 days
	Average walk-away time	40 minutes
	Test selection and arrangement	Manual
	Results per hour*	36
<b>SAMPLING</b>	Sample probe type	Fixed**
	On-board dilution	No
<b>MAINTENANCE</b>	Start-up time	5 minutes
	Daily maintenance time	5 minutes
	Weekly maintenance time	5 minutes
	Monthly maintenance time	10 minutes
<b>SOFTWARE</b>	Software	Built-in
	Operational interface	Touch screen
	Result transfer	RS232C
	Internal storage	800 results
<b>DIMENSIONS</b>	Length	40 cm
	Depth	40 cm
	Height	52 cm
	Weight	29 kg

\* with continuous loading  
 \*\* Cross contamination <0.0001%

AIA-Pack<sup>®</sup>



**TOSOH EUROPE N.V.**  
 sales-marketing.tsen@tosoh.com