



Automated Immunoassay Analysers

AIA[®]-360

The Special One



Automated Immunoassay Analyser AIA[®]-360

Better diagnostics **for better living**

TOSOH BIOSCIENCE

AIA[®]-360

The Special One



METHODOLOGY	Reagent type	AIA-Pack [®]
	Test principle	1 AIA-Pack [®] = 1 test
	Immunochemical determination	Entirely inside AIA-Pack [®]
	Signal detection	Bichromatic Fluorescence Kinetics
CAPACITY	Time to first result	18 minutes
	Sample loading capacity	25
	Number of different tests per sample	4
	Test loading capacity	25 (if 1 per sample)
	Calibration stability	90 days
	Average walk-away time	40 minutes
	Test selection and arrangement	Manual
	Results per hour*	36
SAMPLING	Sample probe type	Fixed**
	On-board dilution	No
MAINTENANCE	Start-up time	5 minutes
	Daily maintenance time	5 minutes
	Weekly maintenance time	5 minutes
	Monthly maintenance time	10 minutes
SOFTWARE	Software	Built-in
	Operational interface	Touch screen
	Result transfer	RS232C
	Internal storage	800 results
DIMENSIONS	Length	40 cm
	Depth	40 cm
	Height	52 cm
	Weight	29 kg

* with continuous loading
 ** Cross contamination <0.0001%

AIA-Pack[®]



TOSOH EUROPE N.V.
 sales-marketing.tsen@tosoh.com